

Plasminogen antibody [9F9-C4]

Cat. No. GTX79688

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human, Bovine

References (1)

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 91 kDa. ([Note](#))

Product Note Recognises Plasminogen/Plasmin

Properties

Form	Liquid
Buffer	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length native protein (purified) (Human).
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX79688 WB Image

WB analysis of various sample lysates using GTX79688 Plasminogen antibody [9F9-C4].

Lane 1 : Plasminogen Non-Reduced

Lane 2 : Plasmin Non-Reduced

Lane 3 : Kringle-I Fusion

Lane 4 : Plasminogen Reduced

Lane 5 : Plasmin Reduced

Lane 6 : MW Markers



For full product information, images and publications, please visit our [website](https://www.genetex.com).