

Plasminogen antibody (FITC)

Cat. No. GTX79698

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	ELISA
Reactivity	Mouse, Rat

Package
500 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Not tested in other applications.	
Calculated MW	91 kDa. (Note)
Product Note	Recognises Plasminogen/Plasmin

PROPERTIES

Form	Liquid
Buffer	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length native protein (purified) (Mouse).
Purification	Protein A purified From serum
Conjugation	Fluorescein isothiocyanate (FITC)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).