

## Factor X antibody (FITC)

Cat. No. GTX79775

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Human

Package  
1 mg

## Applications

## Application Note

Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

## Properties

Form	Liquid
Buffer	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	10 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Affinity purified by Protein G.
Conjugation	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).