

alpha Tubulin antibody [TU-02]

Cat. No. GTX79869

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB
Reactivity	Human, Mouse, Pig

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Product Note This antibody recognizes an epitope on N-terminal structural domain of alpha-tubulin in various species.

Properties

Form	Liquid
Buffer	TBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	microtubule proteins from porcine brain
Purification	Purified by precipitation and chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX79869 WB Image

WB analysis of Jurkat cell lysate (1% laurylmaltoside - non-reduced sample) using GTX79869 alpha Tubulin antibody [TU-02].



For full product information, images and publications, please visit our [website](#).