

LIME antibody [mLIME-05]

Cat. No. GTX79915

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IP
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))**Product Note** This antibody reacts with the cytoplasmic domain of LIME.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein corresponding to amino acids 86-202 of mouse LIME
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX79915 WB Image

WB analysis of mouse thymocyte lysate (1% SDS - non-reduced sample) using GTX79915 LIME antibody [mLIME-05].



For full product information, images and publications, please visit our [website](#).