

Ku80 antibody [MEM-54]

Cat. No. GTX79937

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IP	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	human whole T-lymphocytes
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX79937 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX79937 Ku80 antibody [MEM-54].

Green : Primary antibody

Red : Actin

Blue : DAPI



For full product information, images and publications, please visit our [website](#).