

CD30 antibody [MEM-268] (APC)

Cat. No. GTX79954

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	FCM
Reactivity	Human

Package 100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
FCM	$10~\mu l$ reagent / $100~\mu l$ of whole blood or 10^6 cells in a suspension	
Not tested in other applications.		
Product Note	This antibody recognizes extracellular part of CD30 (Ki-1 antigen).	

Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Expression vector containing CD30 cDNA (booster suspension of THP-1 cell line)
Purification	Purified by size-exclusion chromatography
Conjugation	Allophycocyanin (APC) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 1