

CD99 antibody [MEM-131] (PE)

Cat. No. GTX80016

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Application	FACS
Reactivity	Human

Package
100 test

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension
Not tested in other applications.	

Calculated MW 19 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	TBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	HPB-ALL human peripheral blood leukemia T-cell line
Purification	Purified by size-exclusion chromatography
Conjugation	Phycoerythrin (PE)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX80016 FACS Image

FACS analysis of human peripheral blood using GTX80016 CD99 antibody [MEM-131] (PE).



For full product information, images and publications, please visit our [website](https://www.genetex.com).