

Integrin beta 1 / CD29 antibody [MEM-101A] (APC)

Cat. No. GTX80026

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	FCM
Reactivity	Human, Dog, Pig

Package 100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10 μ l reagent / 100 μ l of whole blood or 10 6 cells in a suspension
Not tested in other applications.	
Product Note	We do not recommend use of this product for Mouse samples.


Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Raji Burkitt's lymphoma cell line
Purification	Purified by size-exclusion chromatography
Conjugation	Allophycocyanin (APC) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 05 Page 1 of 2

DATA IMAGES



GTX80026 FCM Image

FACS analysis of human peripheral blood using GTX80026 Integrin beta 1 / CD29 antibody [MEM-101A] (APC).



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 05 Page 2 of 2