

RN18S1 (18S) Control Primer Pool

Cat. No. GTX80116

Applications	RT-PCR
Reactivity	Human

Package
50 µl

PRODUCT

Summary RT-PCR Primer Pool for RN18S1 (18S)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
RT-PCR	Assay dependent

Not tested in other applications.

Properties

Buffer TE Buffer

Preservative No Preservative

Storage Store at -20°C.

Concentration 100 µM (Please refer to the vial label for the specific concentration.)

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).