

## CPNE7 antibody [CPNE7-01]

Cat. No. GTX80188

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Human      |

References ( 1 )

Package

100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution                                |
|-----------------------------------|-----------------------------------------------------|
| WB                                | Assay dependent                                     |
| Not tested in other applications. |                                                     |
| Calculated MW                     | 70 kDa. ( <a href="#">Note</a> )                    |
| Product Note                      | This antibody recognizes C terminus of human CPNE7. |

## Properties

|               |                                                                                                        |
|---------------|--------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                 |
| Buffer        | PBS                                                                                                    |
| Preservative  | 15mM Sodium azide                                                                                      |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                               |
| Immunogen     | Bacterially produced GST-fused C terminus of human CPNE7.                                              |
| Purification  | Protein A purified                                                                                     |
| Conjugation   | Unconjugated                                                                                           |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX80188 WB Image

WB analysis of nuclear (1) and cytoplasmic fraction (2) of HeLa cell extracts using GTX80188 CPNE7 antibody [CPNE7-01].

Diultion : 1µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).