

E-Cadherin antibody [67A4] (APC)

Cat. No. GTX80259

| Host | Mouse |
|--------------|------------|
| Clonality | Monoclonal |
| Isotype | lgG1 |
| Applications | FCM |
| Reactivity | Human |

Package 100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| FCM $10 \mu l$ reagent / $100 \mu l$ of whole blood or 10^6 cells in a suspension | Suggested dilution | Recommended dilution |
|---|--------------------|---|
| | FCM | 10 μl reagent / 100 μl of whole blood or 10^6 cells in a suspension |

Not tested in other applications.

| Properties | |
|---------------|---|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 15mM Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | T-47D cells |
| Purification | Purified by size-exclusion chromatography |
| Conjugation | Allophycocyanin (APC) Wavelength |
| Note | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
| | Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |



For full product information, images and publications, please visit our website.

Date 2025 / 12 / 12 Page 1 of 1