

## ICAM1 / CD54 antibody [1H4] (FITC)

Cat. No. GTX80290

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	FACS
Reactivity	Human

Reference ( 1 )  
Package  
100 test

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	20 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension
Not tested in other applications.	
Calculated MW	58 kDa. ( <a href="#">Note</a> )
Observed MW	58 or 85~110 kDa. The observed M.W. is based on the following publication. PMID: 24474243

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Raji cells and spleen cells fused with NS1 cells
Purification	Purified by size-exclusion chromatography
Conjugation	Fluorescein isothiocyanate (FITC)


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX80290 FACS Image**

FACS analysis of human peripheral blood using GTX80290 ICAM1 / CD54 antibody [1H4] (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).