

## ICAM1 / CD54 antibody [1H4] (PerCP)

**Cat. No. GTX80292**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	FCM
<b>Reactivity</b>	Human

References ( 1 )

Package

100 test

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension

**Note : Do not need a permeabilization step.**

Not tested in other applications.

**Observed MW (kDa)** 58 or 85~110 kDa. The observed M.W. is based on the following publication. PMID: 24474243

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	15mM Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Raji cells and spleen cells fused with NS1 cells
<b>Purification</b>	Purified by size-exclusion chromatography
<b>Conjugation</b>	Peridinin-chlorophyll proteins (PerCP) <a href="#">Wavelength</a>


**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX80292 FCM Image**

Surface staining of U937 human histiocytic lymphoma cell line with anti-human CD54 (1H4) PE. Total viable cells were used for analysis.



For full product information, images and publications, please visit our [website](#).