

## CD90.2 antibody [1aG4/C5] (PE)

Cat. No. GTX80294-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Application	FACS
Reactivity	Mouse

Reference ( 1 )

Package

100 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	Assay dependent

Not tested in other applications.

Calculated MW 18 kDa. ( [Note](#) )

**Product Note** This antibody recognizes alloantigen Thy-1.2 of mouse Thy-1 (CD90); it does not cross-react with Thy-1.1 alloantigen. We do not recommend use of this product for Human,Rat samples.

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	thymocytes from C3H/Di mice (Thy-1.2 positive)
Purification	Purified by size-exclusion chromatography
Conjugation	Phycoerythrin (PE)

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).