

## PDE4A antibody, C-term

Cat. No. GTX80311

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, IHC
Reactivity	Human, Mouse, Rat

Package

50 µl

## Applications

## Application Note

## Recommended Starting Dilutions:

For WB: Use at 1:2000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	98 kDa. ( <a href="#">Note</a> )
Product Note	The antibody selectively detects proteins corresponding to members of the PDE4A sub-family, in human, mouse & rat

## Properties

Form	Liquid
Buffer	Diluted serum
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Designed from the C-terminal region common to all 4A isoforms
Purification	Unpurified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).