

Cytokeratin 3 antibody [AE5]

Cat. No. GTX80319

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Rabbit, Bovine

Package
50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:1000-1:3000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 64 kDa. ([Note](#))**Product Note** The 64K basic keratin 3 (human, rabbit and bovine)

Properties

Form Liquid**Buffer** PBS**Preservative** 0.1% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Immunogen** Rabbit corneal epithelial keratins**Purification** Protein G purified**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).