

# TRF1 antibody [3H11]

**Cat. No. GTX80331**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Mouse             |
| <b>Clonality</b>    | Monoclonal        |
| <b>Isotype</b>      | IgG1              |
| <b>Applications</b> | WB, ICC/IF, ELISA |
| <b>Reactivity</b>   | Human             |

**Package**

50 µl

## Applications

### Application Note

#### Recommended Starting Dilutions:

For WB: Use at 1:100-1:200

For ICC/IF: Use at 1:100-1:200. Use Hela cells as positive control

Optimal dilutions should be determined experimentally by the researcher.

|                      |                                  |
|----------------------|----------------------------------|
| <b>Calculated MW</b> | 50 kDa. ( <a href="#">Note</a> ) |
| <b>Product Note</b>  | Human TRF1                       |

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.1% Sodium azide                                                                                                                                                                                                          |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | Human TRF1                                                                                                                                                                                                                 |
| <b>Purification</b>  | Protein A                                                                                                                                                                                                                  |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).