

TAP1 antibody [53H8]

Cat. No. GTX80332

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Human

Package

50 µl

Applications

Application Note

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 87 kDa. ([Note](#))

Product Note Human TAP1. Detects a band of approximately 68 kDa

Properties

Form Liquid

Buffer PBS

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Purification Protein A

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).