

Serine Protease Inhibitor 2A (Spi2A) antibody [MoFo29.2]

Cat. No. GTX80339

Host	Armenian Hamster
Clonality	Monoclonal
Isotype	lgG2b
Applications	WB
Reactivity	Mouse

Package $50\,\mu\text{l}$

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
WB	Assay dependent	
Not tested in other applications.		
Product Note	Mouse Spi2A	
Properties		
Form	Liquid	
Buffer	PBS	
Preservative	0.1% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)	
Immunogen	Peptide conjugated to KLH: 406-426 NPERSTNFPNGEGASSQR(C)	
Purification	Protein G	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 03 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com