

PTPRB antibody [122.2]

Cat. No. GTX80343

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB, ICC/IF, IP
Reactivity	Rat

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

Use between 1:1000-1:2500

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	224 kDa. (Note)
Product Note	PTPRB

Properties

Form	Liquid
Buffer	Growth Medium (protein free), concentrated and dialysed against PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Embryonic rat brain proteoglycans
Purification	Concentrated
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).