

Podoplanin antibody [8.1.1]

Cat. No. GTX80352

Host	Golden Syrian Hamster
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IP, IHC
Reactivity	Mouse

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW	17 kDa. (Note)
Product Note	GP38/T1alpha/podoplanin

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Thymic epithelial cell line
Purification	Protein G
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).