

## ORC1 antibody [7A7]

Cat. No. GTX80360

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Application	WB, IP
Reactivity	Human

Package

50 µl

## APPLICATION

## Application Note

## Recommended Starting Dilutions:

For WB: Use at 1:1000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	97 kDa. ( <a href="#">Note</a> )
Product Note	Human ORC1

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ORC1 (hsORC)
Purification	Protein G
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).