

NQO1 antibody [AI80]

Cat. No. GTX80363

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IP
Reactivity	Human, Rat

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:2000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 31 kDa. ([Note](#))

Product Note Human & Rat NQO1. Does not react with mouse NQO1

Properties

Form Liquid

Buffer PBS

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Immunogen Full length purified recombinant human NQO1

Purification Protein G purified

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).