

Nidogen 1 antibody [JF3], G1 Domain

Cat. No. GTX80365

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA
Reactivity	Human, Mouse

Package

50 µl

Applications

Application Note

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 136 kDa. ([Note](#))**Product Note** Recognises native and denatured mouse and human nidogen 1, and blocks binding of nidogen to laminin

Properties

Form	Liquid
Buffer	500mM Acetic acid, neutralized with 2M Trizma
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment corresponding to the G1/Link region of mouse nidogen 1.
Purification	Protein G
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).