

hnRNP L antibody [4D11]

Cat. No. GTX80386

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Human, Mouse, Newt

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:5000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 64 kDa. ([Note](#))

Product Note Human hnRNP L. Detects a band of approximately 65 kDa. It cross reacts with mouse & newt

Properties

Form Liquid

Buffer PBS

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Immunogen Human hnRNP L proteins (from HeLa cells) purified by affinity chromatography on ssDNA agarose

Purification Protein A

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).