

SETDB1 antibody [5H6A12]

Cat. No. GTX80408

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse, Monkey

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 143 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human SETDB1 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx80408 ICC/IF Image

ICC/IF analysis of LOVO cells using GTx80408 SETDB1 antibody [5H6A12].

Green : SETDB1

Red: Actin filaments



GTx80408 WB Image

WB analysis of MCF-7 (1), T47D (2), HEK293 (3), JURKAT (4), NIH3T3 (5), F9 (6), RAW246.7 (7) and Cos7 (8) cell lysate using GTx80408 SETDB1 antibody [5H6A12].



For full product information, images and publications, please visit our [website](https://www.genetex.com).