

BOB1 antibody, N-term

Cat. No. GTX80565

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	1:10-1:50
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 1-27 amino acids from the N-terminal region of human POU2AF1.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx80565 WB Image


WB analysis of Ramos whole cell lysate using GTx80565 BOB1 antibody, N-term.

Loading : 20 µg per lane

Dilution : 1:2000


GTx80565 IHC-P Image

IHC-P analysis of human tonsil using GTx80565 BOB1 antibody, N-term.


GTx80565 FCM Image

FACS analysis of Ramos cells using GTx80565 BOB1 antibody, N-term.

Top histogram : negative control

Bottom histogram : Ramos cells



For full product information, images and publications, please visit our [website](http://www.genetex.com).