

## LMX1B antibody, Internal

Cat. No. GTX80930

|              |              |
|--------------|--------------|
| Host         | Rabbit       |
| Clonality    | Polyclonal   |
| Isotype      | IgG          |
| Applications | WB, IHC-P    |
| Reactivity   | Human, Mouse |

Package  
400 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| IHC-P              | 1:50-1:100           |

Not tested in other applications.

Calculated MW 45 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH conjugated synthetic peptide between 175-202 amino acids from the Central region of human LMX1B.                                                                                                                       |
| Purification  | Protein A purified, followed by peptide affinity purification.                                                                                                                                                             |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx80930 WB Image

WB analysis of mouse heart tissue lysate (35ug/lane) using GTx80930 LMX1B antibody, Internal.



### GTx80930 WB Image

WB analysis of 293 cell lysate (2 ug/lane) either nontransfected (Lane 1) or transiently transfected with the LMX1B (Lane 2) using GTx80930 LMX1B antibody, Internal.



### GTx80930 IHC-P Image

IHC-P analysis of human testis tissue using GTx80930 LMX1B antibody, Internal.



For full product information, images and publications, please visit our [website](https://www.genetex.com).