

## MGST2 antibody, C-term

Cat. No. GTX80967

|              |                |
|--------------|----------------|
| Host         | Rabbit         |
| Clonality    | Polyclonal     |
| Isotype      | IgG            |
| Applications | WB, IHC-P, FCM |
| Reactivity   | Human          |

Package  
400 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| IHC-P              | 1:50-1:100           |
| FCM                | 1:10-1:50            |

Not tested in other applications.

Calculated MW 17 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH conjugated synthetic peptide between 117-147 amino acids from the C-terminal region of human MGST2.                                                                                                                    |
| Purification  | Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.                                                                                                                                           |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX80967 IHC-P Image**


IHC-P analysis of human prostate carcinoma using GTX80967 MGST2 antibody, C-term.


**GTX80967 FCM Image**

FACS analysis of HepG2 cells using GTX80967 MGST2 antibody, C-term.

Green : primary antibody

Blue : negative control


**GTX80967 WB Image**

WB analysis of HepG2 cell lysate (35ug/lane) using GTX80967 MGST2 antibody, C-term.



For full product information, images and publications, please visit our [website](http://www.genetex.com).