

OLR1 / LOX-1 antibody, Internal

Cat. No. GTX80987

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:10-1:50

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 64-92 amino acids from the Central region of human OLR1.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX80987 WB Image


WB analysis of HeLa cell , human liver, and placenta tissue lysate using GTX80987 OLR1 / LOX-1 antibody, Internal.

Loading : 35ug per lane

Dilution : 1:1000


GTX80987 WB Image

WB analysis of 293 cell lysate (2 ug/lane) either nontransfected (Lane 1) or transiently transfected with the OLR1 (Lane 2) using GTX80987 OLR1 / LOX-1 antibody, Internal.


GTX80987 WB Image

WB analysis of HL-60 cell lysate (35ug/lane) using GTX80987 OLR1 / LOX-1 antibody, Internal.



For full product information, images and publications, please visit our [website](http://www.genetex.com).