

ITGBL1 antibody, Internal

Cat. No. GTX81030

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 334-360 amino acids from the Central region of human ITGBL1.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX81030 IHC-P Image


IHC-P analysis of human brain tissue using GTX81030 ITGBL1 antibody, Internal.


GTX81030 FCM Image

FACS analysis of NCI-H460 cells using GTX81030 ITGBL1 antibody, Internal.

Green : primary antibody

Blue : negative control


GTX81030 WB Image

WB analysis of HeLa cell lysate (35ug/lane) using GTX81030 ITGBL1 antibody, Internal.



For full product information, images and publications, please visit our [website](https://www.genetex.com).