

CSGALNACT2 antibody, Internal

Cat. No. GTX81034

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM, IHC
Reactivity	Human

References (2)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100
FCM	1:10-1:50
IHC	Assay dependent

Not tested in other applications.

Calculated MW 63 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 200-229 amino acids from the Central region of human CSGALNACT2.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx81034 WB Image

WB analysis of 293 whole cell lysate using GTx81034 CSGALNACT2 antibody, Internal.

Loading : 20 µg per lane


Dilution : 1:1000


GTx81034 FCM Image


FACS analysis of 293 cells using GTx81034 CSGALNACT2 antibody, Internal.

Green : primary antibody

Blue : negative control


GTx81034 IHC-P Image

IHC-P analysis of human hepatocarcinoma using GTx81034 CSGALNACT2 antibody, Internal.


GTx81034 IHC Image

The data was published in the journal Oncotarget in 2017. [PMID: 28514734](https://pubmed.ncbi.nlm.nih.gov/28514734/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).