

## PXR antibody, Internal

Cat. No. GTX81102

|              |                        |
|--------------|------------------------|
| Host         | Rabbit                 |
| Clonality    | Polyclonal             |
| Isotype      | IgG                    |
| Applications | WB, ICC/IF, IHC-P, FCM |
| Reactivity   | Human                  |

References ( 3 )

Package

400 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| ICC/IF             | 1:10-1:50            |
| IHC-P              | 1:50-1:100           |
| FCM                | 1:10-1:50            |

Not tested in other applications.

Calculated MW 50 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH conjugated synthetic peptide between 100-127 amino acids from the Central region of human NR112.                                                                                                                       |
| Purification  | Protein A purified, followed by peptide affinity purification.                                                                                                                                                             |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX81102 WB Image**


WB analysis of HeLa whole cell lysate using GTX81102 PXR antibody, Internal.

Loading : 20 µg per lane

Dilution : 1:1000


**GTX81102 WB Image**

WB analysis of HL-60 cell lysate (35ug/lane) using GTX81102 PXR antibody, Internal.


**GTX81102 IHC-P Image**

IHC-P analysis of human hepatocarcinoma using GTX81102 PXR antibody, Internal.



For full product information, images and publications, please visit our [website](https://www.genetex.com).