

CCHCR1 antibody, Internal

Cat. No. GTX81194

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:10-1:50
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 89 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 599-627 amino acids from the Central region of human CCHCR1.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx81194 ICC/IF Image

ICC/IF analysis of 293 cells using GTx81194 CCHCR1 antibody, Internal.

Green : CCHCR1


Blue : DAPI


GTx81194 FCM Image

FACS analysis of WiDr cells using GTx81194 CCHCR1 antibody, Internal.

Top histogram : negative control

Bottom histogram : WiDr cells


GTx81194 WB Image

WB analysis of human heart lysate using GTx81194 CCHCR1 antibody, Internal.

Loading : 20 µg per lane

Dilution : 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).