

RUNX2 antibody

Cat. No. GTX81326

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

References (3)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:10-1:50
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 445-474 amino acids of human RUNX2.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX81326 IHC-P Image

IHC-P analysis of human tonsil using GTX81326 RUNX2 antibody.




GTX81326 FCM Image

FACS analysis of NCI-H460 cells using GTX81326 RUNX2 antibody.

Green : primary antibody

Blue : negative control



GTX81326 WB Image

WB analysis of CEM cell lysate (35ug/lane) using GTX81326 RUNX2 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).