

Vitronectin antibody, N-term

Cat. No. GTX81434

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:8000
IHC-P	1:10-1:50
FCM	1:25

Note : Permeabilization : 90% methanol for 10 min

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 65-93 amino acids from the N-terminal region of human VTN.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX81434 WB Image

WB analysis of human plasma using GTX81434 Vitronectin antibody, N-term.

Loading : 20 µg per lane

Dilution : 1:32000




GTX81434 WB Image

WB analysis of human liver lysate using GTX81434 Vitronectin antibody, N-term.

Loading : 20 µg per lane

Dilution : 1:8000



GTX81434 IHC-P Image

IHC-P analysis of human lung carcinoma using GTX81434 Vitronectin antibody, N-term.



For full product information, images and publications, please visit our [website](https://www.genetex.com).