

## HSP47 antibody, C-term

Cat. No. GTX81503

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package  
400 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 46 kDa. ( [Note](#) )

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 390-418 amino acids from the C-terminal region of human SERPINH1.
Purification	Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Conjugation	Unconjugated


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX81503 FCM Image**


FACS analysis of HepG2 cells using GTX81503 HSP47 antibody, C-term.

Green : primary antibody

Blue : negative control


**GTX81503 WB Image**

WB analysis of HeLa cell lysate (35ug/lane) using GTX81503 HSP47 antibody, C-term.


**GTX81503 IHC-P Image**

IHC-P analysis of human testis tissue using GTX81503 HSP47 antibody, C-term.



For full product information, images and publications, please visit our [website](https://www.genetex.com).