

## C9orf156 antibody, C-term

**Cat. No. GTX81642**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human      |

**Package**  
400 µl

### Applications

#### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| IHC-P              | 1:50-1:100           |

Not tested in other applications.

**Calculated MW** 49 kDa. ( [Note](#) )

### Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | KLH conjugated synthetic peptide between 365-391 amino acids from the C-terminal region of human C9orf156.                                                                                                                 |
| <b>Purification</b>  | Protein A purified, followed by peptide affinity purification.                                                                                                                                                             |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

#### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX81642 IHC-P Image

IHC-P analysis of human skin carcinoma using GTX81642 C9orf156 antibody, C-term.



### GTX81642 WB Image

WB analysis of HeLa, CEM, T47D cell lysate (35ug/lane) using GTX81642 C9orf156 antibody, C-term.



For full product information, images and publications, please visit our [website](#).