

ESPN antibody, N-term

Cat. No. GTX81675

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

References (1)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	Assay dependent
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 92 kDa. ([Note](#))

Product Note KO/KD validation is based on published data (PMID: 30618491).

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 17-45 amino acids from the N-terminal region of human ESPN.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx81675 IHC-P Image


IHC-P analysis of human testis tissue using GTx81675 ESPN antibody, N-term.


GTx81675 FCM Image

FACS analysis of 293 cells using GTx81675 ESPN antibody, N-term.

Green : primary antibody

Blue : negative control


GTx81675 WB Image

WB analysis of SK-BR-3 whole cell lysate using GTx81675 ESPN antibody, N-term.

Dilution : 1:1000



For full product information, images and publications, please visit our [website](#).