

Fibrinogen alpha antibody, N-term

Cat. No. GTX81691

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

References (1)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:10-1:50
IHC-P	1:10-1:50

Not tested in other applications.

Calculated MW 95 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 116-144 amino acids from the N-terminal region of human FGA.
Purification	Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX81691 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX81691 Fibrinogen alpha antibody, N-term.


Green : Fibrinogen alpha

Blue : DAPI



GTX81691 WB Image

WB analysis of human blood plasma(35ug/lane) using GTX81691 Fibrinogen alpha antibody, N-term.



GTX81691 IHC-P Image

IHC-P analysis of human breast carcinoma using GTX81691 Fibrinogen alpha antibody, N-term.



For full product information, images and publications, please visit our [website](https://www.genetex.com).