

GARP antibody, Internal

Cat. No. GTX81794

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (1)

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	Assay dependent
IHC-P	1:10-1:50

Not tested in other applications.


Calculated MW 72 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 234-260 amino acids from the Central region of human GARP.
Purification	Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX81794 IHC-P Image

IHC-P analysis of human cerebellum tissue using GTX81794 GARP antibody, Internal.


GTX81794 WB Image

WB analysis of mouse cerebellum tissue lysate (35ug/lane) using GTX81794 GARP antibody, Internal.


GTX81794 WB Image

WB analysis of CEM cell lysate (35ug/lane) using GTX81794 GARP antibody, Internal.



For full product information, images and publications, please visit our [website](https://www.genetex.com).