

FNIP1 antibody

Cat. No. GTX82555

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:1000-2000 in 1% skim milk/TBST after blocking in 5% skim milk/TBST
Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	131 kDa. (Note)
Product Note	Human and mouse FNIP1

Properties

Form	Liquid
Buffer	Serum
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A peptide at the N-terminus of FNIP1, MAPTLFQKLFSKRTGLGAC-amide
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).