

DDDDK tag antibody [FG4R]

Cat. No. GTX82562

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, ELISA, IHC
Reactivity	Species independent

References (4)

★★★★★ Review (1)

Package

100 µg

Applications

Application Note

Recommended Starting Dilutions:

WB (with ECL): Use at a dilution of 1:1000-3000

Immunostaining: Use at a dilution of 1:500 - 1:2000

For best results with other assays (e.g.: Dot, ELISA, IP, etc), please determine optimal working dilution by titration test

Product Note Recognizes the N-terminal, C-terminal or internal DYKDDDDK-tagged fusion proteins

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	DYKDDDDK (FLAG) synthetic peptide conjugated to KLH
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx82562 WB Image

1:1000 (1µg/mL) Ab dilution probed against HEK293 cells transfected with DYKDDDDK-tagged protein vector; untransfected (1) and transfected (2)



For full product information, images and publications, please visit our [website](https://www.genetex.com).