

V5 tag antibody [E10 / V4RR]

Cat. No. GTX82563

| | |
|---------------------|----------------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | WB, ICC/IF, IP, ELISA, IHC |
| Reactivity | Species independent |

 Review (1)

Package
 100 µg

Applications

Application Note**Recommended Starting Dilutions:**

WB (with ECL): Use at a dilution of 1:1000-3000

Immunostaining: Use at a dilution of 1:500 - 1:2000

For best results with other assays (e.g.: Dot, ELISA, IP, etc), please determine optimal working dilution by titration test

Product Note Recognizes V5-tagged fusion proteins

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | No preservatives |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | GKPIPPLLGLDST (V5) synthetic peptide conjugated to KLH |
| Purification | Protein A affinity chromatography from culture supernatant of hybridoma cells grown in a bioreactor |
| Conjugation | Unconjugated |


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


 For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX82563 WB Image**

1:000 (1 μ g/mL) Ab dilution probed against HEK293 cells transfected with V5-tagged protein vector; untransfected (1) and transfected (2)



For full product information, images and publications, please visit our [website](#).