

OLIG1 antibody

Cat. No. GTX82621

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP, IHC
Reactivity	Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:3000
ICC/IF	1:100
IP	1:50
IHC	1:100

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	200mM Tris, 100mM Glycine
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant mouse OLIG1.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX82621 Image



GTX82621 WB Image

WB analysis of neonatal rat brain lysats using GTX82621 OLIG1 antibody.



GTX82621 ICC/IF Image

ICC/IF analysis of primary rat oligodendrocytes using GTX82621 OLIG1 antibody.

Red : Primary antibody

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).