

# Nuclear pore Complex antibody [39C7]

**Cat. No. GTX82625**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Mouse             |
| <b>Clonality</b>    | Monoclonal        |
| <b>Isotype</b>      | IgG1              |
| <b>Applications</b> | ICC/IF            |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µl

## Applications

### Application Note

ICC/IF: 1:50-1:100


## Properties

|                      |                                                                                                                                                                                                                                                                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| <b>Preservative</b>  | 10mM Sodium azide                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>     | Yeast nuclear preparations                                                                                                                                                                                                                                                                                                                  |
| <b>Purification</b>  | Total IgG fraction                                                                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| <b>Note</b>          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX82625 ICC/IF Image



For full product information, images and publications, please visit our [website](#).