

NMDAR2B antibody, N-Term

Cat. No. GTX82630

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC
Reactivity	Rat

Package
100 µl

Applications

Application Note

Suggested starting dilutions are as follows.

WB: 1:1000

IHC: 1:500

IF: 1:500

Not yet tested in other applications. Optimal working dilutions should be determined experimentally by the end user.

Calculated MW 166 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM HEPES, 150mM NaCl, 0.01% BSA, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to amino acid residues from the N-terminal region of the NR2B subunit conjugated to KLH
Purification	Affinity Purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

Anti-NMDA Receptor, NR2B subunit



GTX82630 Image



For full product information, images and publications, please visit our [website](#).