

NMDAR2B antibody

Cat. No. GTX82632

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP, IHC
Reactivity	Human, Mouse, Rat

Package
10 µg

Applications

Application Note

WB: 1:1000

IHC (frozen sections; unpublished observations): 1:1000 to 1:2000

IF: 1:1000 to 1:2000

IP: 3 µl per 200 µg lysate

Calculated MW 166 kDa. ([Note](#))

Properties

Form	Liquid
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Fusion protein from the C-terminal region of the NR2B subunit
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX82632 Image


For full product information, images and publications, please visit our [website](https://www.genetex.com).