

NF-H antibody [NAP4]

Cat. No. GTX82647

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC
Reactivity	Human, Rat

Package
100 µl

Applications

Application Note

Optimal working dilutions should be determined experimentally by the end user.

Calculated MW 112 kDa. ([Note](#))

Properties

Form	Liquid
Preservative	10mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified bovine NF-H
Purification	Total IgG fraction
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES
Anti-Neurofilament H

GTx82647 ICC/IF Image

Immunostaining of human cerebellar cortex showing labeling of NF-H (brown) in basket cell axons surrounding the large Purkinje neurons.

Anti-Neurofilament H

GTx82647 WB Image

Western blot of rat cortex lysate showing specific immunolabeling of the ~200k NF-H protein using Neurofilament NF-H antibody [NAP4] (GTx82647).



For full product information, images and publications, please visit our [website](https://www.genetex.com).